

Schwer Lab

Chemically competent DH5 α bacteria

The day before: Rinse two 1-L Erlenmeyer flasks thoroughly with DI water, then autoclave. Streak *E. coli* DH5 α strain directly from frozen stock onto NON-selective 2X LB agar plate and Carbenicillin or Ampicillin plate (control); incubate streaked plates and one unstreaked NON-selective 2X LB plate (control for contamination) at 37°C for 16 h (3 plates total).

Autoclave microcentrifuge tubes and multistep-pipette tips (for 200 μ l). Prepare TfbI and TfbII. Store solutions and tubes in 4°C cold room O/N.

0.1 M MOPS, pH 7.0

2.09 g MOPS in 80 mL water, pH with NaOH to pH 7.0 (DO NOT overshoot pH!).

Make volume up to 100 mL.

TfBI [30mM CH₃COOK (K-acetate, pH 7), 50 mM MnCl₂, 100 mM KCl, 10 mM CaCl₂, 15% glycerol]

1.5 mL 1 M KOAc

2.5 mL 1 M MnCl₂

5 mL 1 M KCl

0.5 mL 1 M CaCl₂

7.5 g glycerol

pH to 7.0 with glacial acetic acid

Add water to 50 mL. Filter-sterilize (0.2- μ m) in tissue culture hood. Cool to 4°C O/N.

TfBII [10mM MOPS (pH 7 with NaOH), 10mM KCl, 75mM CaCl₂, 15% glycerol]

5 mL 0.1 M MOPS

3.75 mL 1 M CaCl₂

0.5 mL 1 M KCl

7.5 g glycerol

pH to 7.0 with NaOH

Add water to 50 mL. Filter-sterilize in tissue culture hood. Cool to 4°C O/N.

TYM

20 g Tryptone

5 g Yeast Extract

6.8 g NaCl

Bring to 1 L with DI water, autoclave. Add 5 mL of filter-sterilized 1 M MgSO₄ to each 500 mL when solution is below 50°C.

The day of: place six 50-mL conical tubes, sterile multistep pipette tips, and serological pipets (six 25-mL, three 10-mL) on ice.

1. Transfer five well-isolated colonies to 1 mL of TYM (unopened bottle) in sterile 15-mL tube. Vortex at moderate speed; pipette 1 mL culture into 100 mL of TYM in each of two sterile 1-L flasks.
2. Grow cells for ~3.5 h at 37°C (220-250 RPM). OD₆₀₀ should be ~0.6 (1 OD₆₀₀ = ~1 billion cells/mL).
3. Near open flame: transfer cells to sterile, ice-cold 50-mL polypropylene tubes. Cool cultures in ice-water slurry for 10 min.
4. Centrifuge at 2700 x g for 10 min at 4°C (swing-out benchtop centrifuge set to 2°C) to pellet cells.
5. Decant the medium; invert tubes for 1 min to drain medium onto paper towels.

6. Resuspend pellets in 20 mL (for each 100-mL culture) of ice-cold TfBI by swirling and VERY GENTLE tapping (like slowly dragging across a tube rack); can also pipette GENTLY with 25-mL pipette; Pool into one tube.
7. Incubate in ice-water slurry for 15-25 min.
8. Recover cells by centrifugation at 2700 x g for 10 min at 4°C (MUST be 4°C).
9. Pour off buffer; invert tubes for 1 min to drain medium from cell pellet.
10. Resuspend pellet by swirling and VERY GENTLY pipetting in 8 mL of ice-cold TfBI.
11. Incubate in ice-water slurry for 5 min.
12. Aliquot 200 µL into ice-cold sterile microcentrifuge tubes that are placed into cardboard box with holes in bottom on wet ice.
13. After aliquoting all, immerse box in dry ice with 95% ethanol for snap freezing. Remove box and place on dry ice to get rid of excess ethanol. Store box at -80°C.