

## Schwer Lab

### Some buffers for DNA work:

**5X TBE** (89 mM Tris, 89 mM boric acid, 2 mM EDTA)

To 6 Liters of distilled water add:

Tris Base (121.14 g/mol)	432 g
Boric Acid (61.83 g/mol)	220 g
EDTA dihydrate (372.24 g/mol)	29.8 g

Bring total volume to 8 L. Add stir bar and stir o/n at RT.

### 20x E buffer for Southern blotting

10 Liters:

970 g	Tris base (MW 121.1)
1360 g	Sodium acetate trihydrate (MW 136.1)
74 g	EDTA dihydrate (MW 372.2)
500 mL	HCl concentrated (MW 36.5)

Adjust pH to 7.8 - 9.0. Store at RT.

### 20x SSC for Southern blotting

NaCl	1753.2 g
Tri-Na Citrate ( $\text{Na}_3\text{C}_6\text{H}_5\text{O}_7 \cdot 2\text{H}_2\text{O}$ ; 294.1 g/mol)	882.3 g

distilled water to 10 L total volume

### 20x SSCPE

10 Liters:

1400 g	NaCl (MW 58.4)
881.5 g	Sodium citrate dihydrate (MW 294.1)
354 g	Monobasic potassium phosphate (MW 136.09)
74 g	EDTA (MW 372.2)

Adjust pH to 7.2 with NaOH. Can be autoclaved (20 min). Store at RT.

### Standard TE Buffer [10 mM Tris-Cl, pH 8, 1 mM EDTA]

0.5 mL 1 M Tris, pH 8  
100  $\mu\text{L}$  0.5 M EDTA, pH 8  
MilliQ water to 50 mL

### TE Buffer (low EDTA) for genomic DNA [10 mM Tris-Cl, pH 8, 0.1 mM EDTA]

0.5 mL 1 M Tris, pH 8  
10  $\mu\text{L}$  0.5 M EDTA, pH 8  
MilliQ water to 50 mL