

## Schwer Lab

### 5× Direct-loading PCR buffer

Optimized by Annika Carlson

#### REAGENTS

**1 M KCl** (MW =74.55 g/mol)

Dissolve 3.728 g KCl in 30 mL in 50-mL conical tube. Mix until dissolved. Bring volume to 50 mL with ddH<sub>2</sub>O. Store at RT.

**1 M Tris-Cl, pH 9.0** (MW=121.14 g/mol)

Dissolve 6.057 g Tris base into 40 mL ddH<sub>2</sub>O. Adjust the pH to 9.0 by adding concentrated HCl (~2.1-3.0 mL) while stirring at RT. Bring the volume to 50 mL by adding ddH<sub>2</sub>O. Store at RT.

**10% (w/v) IGEPAL CA-630** in ddH<sub>2</sub>O

**1 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>** (MW = 132.14 g/mol)

Dissolve 6.607 g (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub> in 35 mL ddH<sub>2</sub>O. Bring the volume to 50 mL by adding ddH<sub>2</sub>O. Store at RT.

**1 M MgCl<sub>2</sub>•6H<sub>2</sub>O** (MW = 203.30 g/mol)

Dissolve 10.165 g MgCl<sub>2</sub>•6H<sub>2</sub>O in 35 mL ddH<sub>2</sub>O. Bring the volume to 50 mL by adding ddH<sub>2</sub>O. Store at RT.

**Cresol red (C<sub>21</sub>H<sub>18</sub>O<sub>5</sub>S; Sigma 114472-5G; MW = 382.43 g/mol)**

Migrates at about 50-bp DNA fragment size.

**Sucrose: S0389-500G (Sigma)**

#### PROCEDURE

**5× Direct-loading PCR buffer** [30% (w/v) sucrose, 15 mM MgCl<sub>2</sub>, 50 mM KCl, 40 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 0.03% (w/v) Cresol red, 0.25% (w/v) IGEPAL CA-630, 50 mM Tris-Cl, pH 9.0]

Dissolve 3 g sucrose in 5 mL ddH<sub>2</sub>O

add 150 μL 1 M MgCl<sub>2</sub>•6H<sub>2</sub>O

add 500 μL 1 M KCl

add 400 μL 1 M (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>

add 500 μL 1 M Tris-Cl, pH 9.0

add 0.003 g Cresol red

add 250 μL 10% IGEPAL CA-630

Mix well. Centrifuge briefly. Bring total volume to 10 mL with ddH<sub>2</sub>O. Sterilize with 0.22 μm filter. Prepare 1-mL aliquots and store at -20°C.